



Available through 1000Bulbs.com

Hearing Aid Battery LIFE CHART**								
Life (hrs.) = mAh / mA								
	675 640 mAh		13 310 mAh		312 180 mAh		10 105 mAh	
Battery Current (ma)	Hours	Days*	Hours	Days*	Hours	Days*	Hours	Days*
0.2							525	32.8
0.3					600	37.5	350	21.9
0.4					450	28.1	263	16.4
0.5			620	38.8	360	22.5	210	13.1
0.6			517	32.3	300	18.8	175	10.9
0.7			443	27.7	257	16.1	150	9.4
0.8			388	24.2	225	14.1	131	8.2
0.9	711	44.4	344	21.5	200	12.5	117	7.3
1.0	640	40.0	310	19.4	180	11.3	105	6.6
1.2	533	33.3	258	16.1	150	9.4	88	5.5
1.4	457	28.6	221	13.8	129	8.0	75	4.7
1.8	356	22.2	172	10.8	100	6.3	58	3.6
2.0	320	20.0	155	9.7	90	5.6	53	3.3
2.2	291	18.2	141	8.8	82	5.1	48	3.0
2.4	267	16.7	129	8.1	75	4.7	44	2.7
2.6	246	15.4	119	7.5	69	4.3	40	2.5
2.8	229	14.3	111	6.9	64	4.0	38	2.3
3.0	213	13.3	103	6.5	60	3.8	35	2.2

\* Days of service based on: 16 hours use per day  
 \*\* Actual battery life will vary by individual and hearing aid

IMPORTANT NOTE: Maximum battery life with any zinc air battery will not typically exceed 45 days regardless of what is indicated by the chart. The battery will dry out and no longer be effective.